



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone (800) 325-5832 (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

JumpStart™ Taq DNA Polymerase

With 10× Reaction Buffer containing MgCl_2

Product Code **D9307**

Technical Bulletin No. MB-395

TECHNICAL BULLETIN

Product Description

JumpStart™ Taq DNA polymerase (JumpStart Taq) is an optimized blend of Sigma's high-performance Taq DNA Polymerase and JumpStart™ Taq antibody. In JumpStart Taq DNA polymerase, the Taq activity is inactivated by combining the enzyme with JumpStart Taq antibody, a neutralizing monoclonal antibody to Taq DNA polymerase. The antibody inactivation provides a simple, efficient procedure for hot start PCR[†]. Hot start PCR can significantly improve the results of DNA amplifications by reducing the generation of nonspecific amplification products and primer-dimer artifacts. When used in PCR, JumpStart Taq DNA polymerase is inactive at low (room) temperature. When the temperature is raised above 70 °C in the first denaturation step of the cycling process, the complex dissociates and the polymerase becomes fully active.

Typical applications for JumpStart Taq include PCR reactions involving one or more of the following: complex genomic or cDNA templates, very low-copy-number targets, large number of thermal cycles (>35) and multiple primer pairs in the same reaction tube.

The enzyme is supplied at 2.5 units/ μl and comes with an optimized 10×-reaction buffer containing magnesium chloride.

Unit Definition: One unit incorporates 10 nmol of total deoxyribonucleoside-triphosphates into acid precipitable DNA in 30 minutes at 74 °C.

Reagents Provided

- JumpStart Taq DNA polymerase, Product Code D6558
2.5 Units/ μl in 75 mM KCl, 15 mM Tris-HCl, pH 7.5, 50% glycerol, 0.05 mM EDTA, 0.5 mM DTT, 0.25% Tween[®] 20, and 0.25% IGEPAL[®] CA-630. Provided as 50, 250, or 1500 units.
- 10× PCR Buffer, Product Code P 2192
100 mM Tris-HCl (pH 8.3 at 25 °C), 500 mM KCl, 15 mM MgCl_2 , 0.01% (w/v) gelatin. Provided in 1.5 ml vials.

Reagents and Equipment Required But Not Provided
(Sigma Product Numbers have been given where appropriate.)

- 10 mM dATP, Product Code D6920
- 10 mM dCTP, Product Code D7045
- 10 mM dGTP, Product Code D7170
- 10 mM TTP, Product Code T7791
or
- Deoxynucleotide mix, Product Code D7295
containing 10 mM dATP, 10 mM dCTP, 10 mM dGTP, 10 mM TTP
- PCR reagent water, Product Code W1754
- Mineral oil, Product Code M8662
- 0.2 and 0.5 ml PCR microcentrifuge tubes, thin-walled, Product Codes P3114 and P3364
- Thermal cycler
- Primers
- DNA to be amplified
- Chloroform, Product Code C7559 (optional)

Precautions and Disclaimer

JumpStart Taq DNA polymerase is for R&D use only. Not for drug, household or other uses. Refer to Material Safety Data Sheet.

Storage/Stability

Store at -20°C .

JumpStart Taq DNA polymerase may be stored at 4°C for up to 3 months and at -20°C for up to six months. In the supplied storage buffer and at the supplied concentration, JumpStart Taq will not freeze at -20°C . Freezing JumpStart Taq at storage below -20°C is not recommended. Repeated freeze-thaw of JumpStart Taq may adversely affect its function.

Procedure

Note: The use of DMSO or formamide with JumpStart Taq DNA polymerase is not recommended due to interference with the enzyme-antibody complex. Other cosolvents, solutes (e.g., salts) and extremes in pH or other reaction conditions may reduce the affinity of the JumpStart Taq antibody for the Taq polymerase and thereby compromise its effectiveness.

Preparation of PCR Master Mix and Thermal Cycling Parameters

Because the Taq polymerase is a magnesium ion-dependent enzyme, the optimal conditions for the concentration of Taq, template DNA, primers, and MgCl_2 will depend on the system being utilized. Individual component. This is especially true for the JumpStart Taq, cycling parameters, and the MgCl_2 concentration. It is recommended the enzyme and the MgCl_2 be titrated to determine the optimal efficiency.

To minimize tube-to-tube variation, preparation of a PCR master mix with JumpStart Taq is recommended. The amount prepared should be based on the number of PCR reactions to be performed.

- For a single reaction, add the following reagents to a 0.2 or 0.5 ml microcentrifuge tube in the following order:

Amount	Component	Final Concentration
- μl	Water	-
5 μl	10 \times PCR Buffer	1 \times
1 μl^*	10 mM dATP	200 μM
1 μl^*	10 mM dCTP	200 μM
1 μl^*	10 mM dGTP	200 μM
1 μl^*	10 mM TTP	200 μM
- μl	Primers	0.1-0.5 μM
1 μl	JumpStart Taq DNA Polymerase	0.05 units/ μl
- μl	Template DNA (typically 10 ng)	200 pg/ μl
50 μl	Total reaction volume	

*The individual nucleotides (1 μl of each 10 mM solution, 4 μl total) may be replaced by 1 μl of Deoxynucleotide Mix, Product Code D 7295.

- Mix gently and briefly centrifuge to collect all solution at the bottom of the tube.
- Add 50 μl of mineral oil to the top of each tube to prevent evaporation (optional, depending on the model of thermal cycler).
- Amplification parameters will vary depending on the primers and the thermal cycler used. It may be necessary to optimize the system for individual primers, template, and thermal cycler.

Typical cycling parameters:

Initial denaturation 94°C 1 min

25-35 cycles:

Denaturation	94°C	30 sec
Annealing	55°C to 68°C	30 sec
Extension	72°C	1 min (minimum)*
Final extension:	72°C	1 min (minimum)*
Hold	4°C	

* 1 minute minimum or 1 minute per kb expected amplicon.

5. The amplified DNA can be evaluated by agarose gel electrophoresis and subsequent ethidium bromide staining. Mineral oil overlay may be removed by a single chloroform extraction (1:1), recovering the aqueous phase.

JumpStart Taq DNA Polymerase Troubleshooting Guide

Problem	Suggestion
No reduction of nonspecific products is observed when using JumpStart Taq	<p>Test the PCR system using a manual hot start method.</p> <p>The use of DMSO or formamide with JumpStart Taq is not recommended due to interference with the enzyme-antibody complex. Other cosolvents, solutes (e.g., salts) and extremes in pH or other reaction conditions may reduce the affinity of the JumpStart Taq antibody for the Taq polymerase and thereby compromise its effectiveness.</p>
Both the JumpStart Taq PCR and the manual hot start PCR yield multiple nonspecific products	<p>Raise the annealing temperature in 2-3 °C increments. Raising the temperature improves the specificity of binding by the primers, however, it may result in reduced binding and extension of the primers. If raising the annealing temperature causes a reduced yield of the specific product with only a proportional reduction of side reaction products, it may be necessary to redesign the primers.¹</p> <p>Take special precautions to avoid crossover contamination of PCR reactions with both specific and nonspecific PCR products, including primer-dimer artifacts.²</p>
The JumpStart Taq PCR yields more nonspecific products than manual hot start PCR	<p>Titration of JumpStart Taq may be necessary to optimize the PCR reaction conditions, especially if the conditions vary from those described in this document. In this case, start with a working solution that has a two- to four-fold higher concentration of JumpStart Taq than recommended.</p>
The yield of specific product is low using JumpStart Taq.	<p>Increase the reaction volume to 150 µl or more.</p> <p>Increase the number of amplification cycles. If currently using 25-30 cycles, increase the cycle number to 35-40. This should increase yields without significantly increasing side reaction products.</p> <p>Modify the reaction conditions and/or selection of PCR targets to obtain greater opportunities for PCR priming. For example, increase the denaturation time up to 1-1.5 minutes and/or increase the denaturation temperature to as high as 95 °C to overcome denaturation difficulties.</p> <p>The use of DMSO or formamide with JumpStart Taq is not recommended due to interference with the enzyme-antibody complex.</p>

References

1. Huang, L. M. and Jeang, K.-T. *BioTechniques* **16**:242-246 (1994)
2. Kwok, S. and Higuchi, R. *Nature* **339**:237-238 (1989)

General References

- Griffin, H. G. and Griffin, A. M. (Eds.) *PCR Technology: Current Innovations*, CRC Press, 1994. (Product Code Z35,749-9)
- Innis, M. A., et al. (Eds.), *PCR Strategies*, Academic Press, New York (1995). (Product Code Z36,445-2)
- Innis, M., et al. (Eds.), *PCR Protocols: A Guide to Methods and Applications*, Academic Press, San Diego, California (1990). (Product Code P 8177)
- Newton, C. R. (Ed.), *PCR: Essential Data*, John Wiley & Sons, New York (1995). (Product Code Z36,491-6)
- Sambrook, J. F., et al., *Molecular Cloning: A Laboratory Manual*, Third Edition, Cold Spring Harbor Laboratory Press, New York (2000). (Product Code M 8265)

†The PCR process is covered by patents owned by Hoffman-LaRoche, Inc.

NOTICE TO PURCHASER: LIMITED LICENSE

A license under US Patents 4,683,202, 4,683,195, 4,965,188, and 5,075,216 or their foreign counterparts, owned by Roche Molecular Systems and F. Hoffmann-La Roche Ltd ("Roche"), has an up-front fee component and a running-royalty component. The purchase price of this product includes limited, nontransferable rights under the running-royalty component to use only this amount of the product to practice the Polymerase Chain Reaction ("PCR") and related processes described in said patents solely for the research and development activities of the purchaser when this product is used in conjunction with a thermal cycler whose use is covered by the up-front fee component. Rights to the up-front fee component must be obtained by the end user in order to have a complete license to use this product in the PCR process. These rights under the up-front fee component may be purchased from Applied Biosystems or obtained by purchasing an Authorized Thermal Cycler. No right to perform or offer commercial services of any kind using PCR, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is hereby granted by implication or estoppel. Further information on purchasing licenses to practice the PCR Process may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or the Licensing Department at Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.

JumpStart and JumpStart Taq Antibody are licensed under U.S. Patent No. 5,338,671 and 5,587,287 and corresponding patents in other countries.

JumpStart™ Is a trademark of Sigma-Aldrich Biotechnology LP.
Tween® is a trade mark of Uniqema, a business unit of ICI Americas Inc.
IGEPAL® is a registered trade mark of Rhone-Poulenc LP

6/2002

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.